Liquid Antithrombin - 0020002500

Intended use

Automated chromogenic assay for the quantitative determination of Antithrombin in human citrated plasma on the IL Coagulation Systems.

Summary and principle

Antithrombin (AT) or Heparin Cofactor I is the major inhibitor of blood coagulation and is essential for effective heparin therapy. By inhibiting the coagulation proteases, especially thrombin, FXa and FIXa, AT prevents uncontrolled coagulation and thrombosis. Antithrombin deficiency is associated with a high risk of thromboembolic disorders.^{1,2,3}

Liquid Antithrombin can be used to exclude or diagnose hereditary deficiency^{4,5} in patients with a tendency toward thromboembolism, in pre-operative stages, before prescription of oral contraceptives, DIC⁶, nephrotic syndromes, liver diseases⁷ and in therapy with heparin or antithrombin concentrates.^{8,9} The Liquid Antithrombin kit is as an assay based on a synthetic chromogenic substrate and on FXa inactivation.^{10,11} As a consequence, the method is specific and not influenced by Heparin Cofactor II. Antithrombin levels in patient plasma are measured automatically on IL Coagulation Systems in two stages:

1. Incubation of the plasma with the Factor Xa reagent in the presence of an excess of heparin.
2. Quantification of the residual FXa activity with a synthetic chromogenic substrate. The paranitroaniline released is monitored kinetically at 405 nm and is inversely proportional to the Antithrombin level in the test sample.

Composition

The Liquid Antithrombin kit consists of:

- [S] Chromogenic substrate** (Cat. No. 0020002520): 2 x 2 mL vials of the chromogenic substrate S-2765, N- α -Z-D-Arg-Gly-Arg-pNA2HCl (6 mg/vial), surfactant and buffer.
- [E] Factor Xa reagent** (Cat. No. 0020002510): 4 x 4 mL vials of a solution containing bovine Factor Xa (40 nkat/vial), heparin, buffer, bovine serum albumin and preservatives.

PRECAUTIONS AND WARNINGS:

Factor Xa reagent:

- Danger**
- Hazard class:** Resp. Sens. 1, H334
- Hazard statements:** H334: May cause allergy or asthma symptoms or breathing difficulties if inhaled.

Precautionary statements: P261: Avoid breathing vapors/spray. P284: [In case of inadequate ventilation] wear respiratory protection. P304 + P340: IF INHALED: Remove person to fresh air and keep comfortable for breathing. P342 + P311: If experiencing respiratory symptoms: Call a POISON CENTER/doctor. P501: Dispose of contents/container in accordance with local/regional/national/international regulation.

Supplemental hazard information: Contains Factor Xa. Up to 3.4% of the mixture consists of component of unknown acute toxicity (oral, dermal, inhalation) for the human health and unknown hazard to the aquatic environment.

Chromogenic Substrate:

Hazard class: none

Hazard statements: none

Precautionary statements: none

Supplemental hazard information: none

Contains bovine material. All donor animals were sourced from BSE-free herds. The cattle received ante- and post-mortem health inspection by a veterinarian, and they were apparently free from infectious and contagious material. However, the material should be treated as potentially infectious. This product is For *in vitro* Diagnostic Use.

Preparation

Chromogenic substrate: Invert to mix before use.

Factor Xa reagent: Invert to mix before use.

Reagent storage and stability

Unopened reagents are stable until the expiration date shown on the vial when stored at 2-8°C.

Chromogenic substrate - Opened reagent is stable: 5 weeks at 2-8°C in the original vial, 48 hours at 15°C on the ACL ELITE®/ELITE PRO/8/9/10000, ACL Futura and ACL Advance Systems or 8 hours at 15°C on the ACL Classic Systems. Do not freeze.

Factor Xa reagent - Opened reagent is stable: 5 weeks at 2-8°C in the original vial, 48 hours at 15°C on the ACL ELITE/ELITE PRO/8/9/10000, ACL Futura and ACL Advance Systems or 8 hours at 15°C on the ACL Classic Systems. Do not freeze.

For optimal stability remove reagents from the system and store them at 2-8°C in the original vial.

Instrument/test procedures

Refer to the IL instrument's Operator's Manual and/or Application Manual for the complete assay procedure instructions.

Specimen collection and preparation

Nine parts of freshly drawn venous blood are collected into one part trisodium citrate. Refer to CLSI Document H21-A5 for further instructions on specimen collection, handling and storage.¹²

Additional reagents and control plasmas

The following are not supplied with the kit and must be purchased separately.

| | Americas and Pacific Rim Cat. No. | Europe Cat. No. |
|------------------------------|--------------------------------------|--------------------|
| Calibration plasma | 0020003700 | 0020003700 |
| Normal Control | 0020003120/0020003110 | 0020003110 |
| Special Test Control Level 1 | 0020011000 | 0020011000 |
| Special Test Control Level 2 | 0020012000 | 0020012000 |
| Low Abnormal Control | 0020003220/0020003210 | 0020003210 |
| Factor diluent | 0009757600 | 0009757600 |
| Cleaning solution | 0009831700 | 0009831700 |
| Cleaning agent | 0009832700 | 0009832700 |

Quality control

Normal and abnormal controls are recommended for a complete quality control program.¹³ Normal Control, Low Abnormal Control, and Special Test Control Level 1 and 2 are designed for this program. Each laboratory should establish its own mean and standard deviation and should establish a quality control program to monitor laboratory testing. Controls should be analyzed at least once every 8 hour shift in accordance with good laboratory practice. Refer to the instrument's Operator's Manual for additional information. Refer to Westgard *et al* for identification and resolution of out-of-control situations.¹⁴

ENGLISH - Insert revision 03/2016

Results

Antithrombin results are reported in activity (%). Refer to the instrument's Operator's Manual for additional information.

Limitations/interfering substances

Antithrombin results on IL Coagulation Systems are not affected by heparin (UF heparin or LMW heparin) up to 4 U/mL, α_1 -Antitrypsin up to 4 mg/mL, α_2 -macroglobulin up to 7 mg/mL, Heparin Cofactor II up to 4 U/mL, hemoglobin up to 150 mg/dL, bilirubin up to 40 mg/dL and triglycerides up to 500 mg/dL.

Expected values

Antithrombin activity levels in healthy individuals are approximately in the range of 83 - 128%. Antithrombin levels are low in neonates/infants and increase to adult levels by approximately 1 year of age; levels are then slightly higher than in adults up to age 16 year.¹⁵ Due to many variables which may affect results, each laboratory should establish its own normal range.

Performance characteristics

Precision:

Within run and total (run to run and day to day) precision was assessed over multiple runs using both normal and abnormal samples.

| ACL Classic (100-7000) | Mean (% activity) | CV % (Within run) | CV % (Total) |
|-------------------------------|-------------------|-------------------|--------------|
| Normal Control | 102,6 | 2,54 | 3,51 |
| Special Test Control Level 1 | 53,6 | 3,27 | 5,47 |
| Special Test Control Level 2 | 22,1 | 7,54 | 11,52 |
| ACL 8/9/10000/ELITE/ELITE PRO | Mean (% activity) | CV % (Within run) | CV % (Total) |
| Normal Control | 114,0 | 1,63 | 2,47 |
| Low Abnormal Control | 36,6 | 2,98 | 3,83 |
| High Abnormal Control | 24,1 | 2,33 | 4,86 |
| ACL Futura/ACL Advance | Mean (% activity) | CV % (Within run) | CV % (Total) |
| Normal Control | 104,6 | 2,10 | 3,88 |
| Special Test Control Level 1 | 54,2 | 2,47 | 3,82 |
| Special Test Control Level 2 | 25,2 | 12,10 | 13,13 |

Correlation:

| System | n | slope | intercept | r | Comparative method |
|------------------------|-----|-------|-----------|-------|----------------------|
| ACL Classic (100-7000) | 118 | 0,99 | -2,35 | 0,995 | HemosIL Antithrombin |
| ACL Futura/Advance | 118 | 1,01 | -2,22 | 0,994 | HemosIL Antithrombin |

In an additional clinical study comparing Liquid Antithrombin to Antithrombin (Cat. No. 0020008900), samples from 64 patients (26 normals, 11 undergoing heparin therapy, 5 undergoing OAT, 4 undergoing antithrombin substitution therapy, 5 with congenital AT deficiency, 5 with liver disease, 4 with DIC and 4 others with high AT) were evaluated. The correlation (r) was 0.963 on the ACL Family and 0.997 on the ACL Futura/ACL Advance. The precision and correlation results were obtained using specific lots of reagents and controls.

Linearity:

| System | ACL Classic (100-7000) and ACL Futura/ACL Advance | 10 -150 (% activity) |
|--------|---|----------------------|
| | | |

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- Wells PS, Blajchman MA, Henderson P The prevalence of Antithrombin Deficiency in Healthy Blood Donors. A Cross-Sectional Study. Am. J. Hematol. 1994; 4: 321-324.
- Tait FC, Walker ID, Perry DJ, Carrell RW, Islam SIA, McCall F, Mitchell R, Davidson JF. Prevalence of Antithrombin III Deficiency Subtypes in 4000 Healthy Blood Donors. Thromb. Haemost. 1991; 65: 534 (abst).
- Harper PL, Luddington RJ, Daly M. The Incidence of Dysfunctional Antithrombin Variants: Four Cases in 210 Patients with Thromboembolic Disease, Br. J. Haematol. 1991; 77: 360-364.
- Perry DJ. Antithrombin and Its Inherited Deficiencies. Blood. 1994; 8: 37-55.
- Pabinger I, Brucker S, Kyrle PA. Hereditary Deficiency of Antithrombin III, Protein C and Protein S. Prevalence in Patients with a History of Venous Thrombosis and Criteria for Rational Patient Screening, Blood Coag. Fibrinolys. 1992; 3: 547-553.
- Fourrier F, Chopin C, Goudemand J, Henryx S, Caron C, Rime A, Marey A, Lestavel P. Septic Shock, Multigang Failure, and Disseminated Intravascular Coagulation. Compared Patterns of Antithrombin III, Protein C and Protein S Deficiencies, Chest. 1993; 101: 816-823.

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